

INTRODUCING THE NEW ON-Q* PUMP WITH SELECT-A-FLOW*

Non-narcotic, postop pain relief just got more personal.



THE ON-Q* PUMP WITH SELECT-A-FLOW* VARIABLE RATE CONTROLLER, designed with physicians' input, is the optimal postop pain management solution. It has the same continuous non-narcotic benefits you have come to rely on from our fixed flow therapy with the added benefits of the following:

- Titratable, which allows for personalization to each patient's pain relief needs from the beginning to the end of their therapy
 - Can be dialed up or down based on patients' need during physical therapy or increased breakthrough pain
- Increases ability for facility protocol standardization
- Lightens pharmacy load by minimizing number of SKUs and associated filling needs



NEWLY ENHANCED DESIGN

Flow rate prominently displayed in window

Improved low profile and ergonomic design

Arrows indicate that knob is between rate settings



Case will not close if dialed between rate settings

| Model | Fill Volume | Flow Rate |
|--------|-------------|--------------------------------|
| SAF01 | 270 ml | 1-7 ml/hr |
| CB004 | 400 ml | 2-14 ml/hr |
| CB4007 | 400 ml | 1-7 ml/hr (dual infusion) |
| CB005 | 400 ml | 1-7 ml/hr + 5 ml/60 min bolus |
| CB006 | 400 ml | 2-14 ml/hr + 5 ml/30 min bolus |
| CB6004 | 600 ml | 2-14 ml/hr |
| CB6007 | 600 ml | 1-7 ml/hr (dual infusion) |

Ask your Avanos rep how you can improve your patients' post op pain recovery protocols with the new **ON-Q* Pump with Select-A-Flow***
1-800-448-3569



There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings** and **Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.avanospainmanagement.com for additional product safety Technical Bulletins.

AVANOS

avanospainmanagement.com



*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
 © 2018 AVNS. All rights reserved.
 MK-00733 02/2015