

# QUIKBLOC\* OVER-THE-NEEDLE CATHETER SETS




## PERFORM CONTINUOUS PERIPHERAL NERVE BLOCKS IN A SINGLE SHOT

A unique catheter over echogenic needle design to perform a **one-step** continuous peripheral nerve block, saving time compared to traditional through-the-needle techniques.

- Minimal catheter dislodgement – 30X more holding power
- Reduced likelihood of leakage
- Custom-fit catheter securement device
- Available in non-stimulating and stimulating configurations



# QUIKBLOC\* SET COMPONENTS



Injection & Extension Tubing Set



ON-Q\* Catheter Securement Device



Dressings & Adhesive Bandages

PRODUCT CODE	NEEDLE				CATHETER		PKG/CASE	NRFit™
	NEEDLE TIP (FACET BEVEL)	GAUGE	USABLE LENGTH (MM/IN)	STIMULATING CABLE	INFUSION SET	GAUGE		
<b>ON-Q* QUIKBLOC* OVER-THE-NEEDLE CATHETER SETS</b>								
QB10016SG	•	20	100/4		•	16	3	5
QB10016SG-N	•	20	100/4		•	16	3	5 •
QB10016SGC	•	20	100/4	•	•	16	3	5
QB10016SGC-N	•	20	100/4	•	•	16	3	5 •
QB15016SG	•	20	150/6		•	16	5	5
QB15016SG-N	•	20	150/6		•	16	5	5 •
QB15016SGC	•	20	150/6	•	•	16	5	5
QB15016SGC-N	•	20	150/6	•	•	16	5	5 •

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings** and **Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to [www.avanospainmanagement.com](http://www.avanospainmanagement.com) for additional product safety Technical Bulletins.

\*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2020 AVNS. All rights reserved. COPY-04450 05/20