


AVANOS SINGLE SHOT NEEDLES

From stimulating to echogenic, Avanos offers a wide range of Single Shot peripheral nerve block needles to serve a variety of physician preferences and needs.

PROBLOC* II SINGLE SHOT NEEDLES

ProBloc* II is our standard single shot needle for peripheral nerve blocks. With a 30° block bevel, this sleeve-coated needle comes in various lengths and non-stimulating and stimulating options to suit physician preferences.



30° block bevel needle tip



Graduated depth markings



Sleeve-coated with exposed tip




Available in stimulating and non-stimulating versions




1.2" (30mm)
3.5" (80mm)
4" (100mm)
6" (150mm)
options

ECHOBRIGHT* ECHOGENIC SINGLE SHOT NEEDLES

Our enhanced single shot needle, the EchoBright* line of products has specially designed, 360° laser etched echogenic dimples to improve needle visibility for ultrasound-guided peripheral nerve blocks.



21° Facet or Tuohy bevel needle tips



Graduated depth markings




Powder-coated with exposed tip



Available in stimulating and non-stimulating versions



2" (50mm)
4" (100mm)
6" (150mm)
options




Unique dimpling design for improved visualization under ultrasound



360° echogenic pattern




NEEDLE TIP CONFIGURATIONS



30° Block Bevel



21° Facet Bevel



Tuohy Bevel

| Product Code | Needle | | | | | | Pkg/Case | NRFit™ |
|--|---|--|-------|-----------------------|-------------------|--------------|----------|--------|
| | Tip B = Block Bevel F = Facet Bevel T = Tuohy | SL = Sleeve Insulated E = Echogenic | Gauge | Usable Length (mm/in) | Stimulating Cable | Infusion Set | | |
| ECHOBRIGHT* ECHOGENIC SINGLE SHOT NEEDLES | | | | | | | | |
| EBL22050SG | F | E | 22 | 50/2 | | • | 10 | |
| EBL22050SG-N | F | E | 22 | 50/2 | | • | 10 | • |
| EBL22050SGC | F | E | 22 | 50/2 | • | • | 10 | |
| EBL22050SGC-N | F | E | 22 | 50/2 | • | • | 10 | • |
| EBL20100SG | F | E | 20 | 100/4 | | • | 10 | |
| EBL20100SG-N | F | E | 20 | 100/4 | | • | 10 | • |
| EBL20100SGC | F | E | 20 | 100/4 | • | • | 10 | |
| EBL20100SGC-N | F | E | 20 | 100/4 | • | • | 10 | • |
| EBL20150SG | F | E | 20 | 150/6 | | • | 10 | |
| EBL20150SGC | F | E | 20 | 150/6 | • | • | 10 | |
| EBL18050TG | T | E | 18 | 50/2 | | • | 10 | |
| EBL18050TG-N | T | E | 18 | 50/2 | | • | 10 | • |
| EBL18050TGC | T | E | 18 | 50/2 | • | • | 10 | |
| EBL18050TGC-N | T | E | 18 | 50/2 | • | • | 10 | • |
| EBL18100TG | T | E | 18 | 100/4 | | • | 10 | |
| EBL18100TG-N | T | E | 18 | 100/4 | | • | 10 | • |
| EBL18100TGC | T | E | 18 | 100/4 | • | • | 10 | |
| EBL18100TGC-N | T | E | 18 | 100/4 | • | • | 10 | • |
| EBL18150TG | T | E | 18 | 150/6 | | • | 10 | |
| EBL18150TG-N | T | E | 18 | 150/6 | | • | 10 | • |
| EBL18150TGC | T | E | 18 | 150/6 | • | • | 10 | |
| EBL18150TGC-N | T | E | 18 | 150/6 | • | • | 10 | • |
| PROBLOC* II SINGLE SHOT NEEDLES | | | | | | | | |
| HN1-40 | B | SL | 22 | 30/1.2 | | | 10 | |
| HN1-80 | B | SL | 22 | 80/3.15 | | | 10 | |
| HN2-40 | B | SL | 22 | 30/1.2 | | • | 10 | |
| HN2-80 | B | SL | 22 | 80/3.15 | | • | 10 | |
| HN3S-40 | B | SL | 22 | 30/1.2 | • | • | 10 | |
| HN3S-40-N | B | SL | 22 | 30/1.2 | • | • | 10 | • |
| HN3S-80 | B | SL | 22 | 80/3.15 | • | • | 10 | |
| HN3S-80-N | B | SL | 22 | 80/3.15 | • | • | 10 | • |
| HN3S-100 | B | SL | 20 | 100/4 | • | • | 10 | |
| HN3S-150 | B | SL | 20 | 150/6 | • | • | 10 | |

For more information for clinicians and product code ordering information, please visit: www.avanosmedicaldevices.com
Call 1-800-448-3569 in the United States and Canada.