

ΔVΔNOS

The **ambIT® Kids PIB Pain Control System** provides a unique approach to meet the needs of pediatric pain management. Parents rely on doctors for their child's safety and physicians rely on ambIT.

- Kid friendly
- ambIT allows parental control
- Non-addictive and opioid-free
- Faster recovery
- Quick set-up; less than 5 minutes



The ambIT® Kids PIB product portfolio is part of a larger ON-Q* offering of products, including a full range of pumps, needles, catheters, trays and additional accessories, along with a comprehensive training, education and support system to provide effective pain relief for your patients.



ambIT® Kids PIB Infusion Therapy System

Volumetric Accuracy	+/- 6%
Dose Volumes	0 to 50 ml in 0.1 ml steps
Dose Interval Times (hr:min)	00:01 to 24:00 in 1 minute increments
Basal Flow Rates	0 to 20 ml/hr in 0.1 ml/hr steps
Stroke Volume	50 microliters
Bolus Delivery Rate†	210 ml/hr
Bolus Dose Volume	0 to 50 ml in 0.1 ml steps
Bolus Lockout Times (hr:min)	00:01 to 24:00 in 1 minute increments
Volumes to be Infused (ml)	0 to 1000
Pump Mechanism	Microprocessor Controlled Rotary-Peristaltic
Dimensions	2.16 in. x 1.4 in. x 6.875 in. (55 mm x 36 mm x 175 mm)
Weight	6.4 ounces (181.4 grams) with batteries
Power Supply	2 AA 1.5V batteries
Alarms & Alerts	<ul style="list-style-type: none"> • Run Indicator Light • Bolus Infusing • Pause Indicator • Occlusion Downstream • Detached Cassette • Lockout Time Remaining • Low Battery • Dead Battery • Malfunction • Boluses Requested • Infusion Complete • Local Bolus Button Off/On
Reports	<ul style="list-style-type: none"> • Volume Infused • Boluses Requested • Boluses Delivered • Elapsed Time
Delivery Profiles	<ul style="list-style-type: none"> • PIB+PCA • PCA • PIB

† Rate varies

ORDERING

844-4-AVANOS
PainManagement@avanos.com

CUSTOMER FEEDBACK

844-4-AVANOS (Option 1, 3)
PIQ@avanos.com

TECHNICAL SUPPORT

800-444-2728
24/7 Nurse Hotline

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings and Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.avanos.com for additional product safety Technical Bulletins.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2020 AVNS. All rights reserved. COPY-04618 07/2020