

TAKE ADVANTAGE OF THE AMBIT* PIB-IL PUMP AND AMBIT*+ PROGRAM

Discover more ways to bring value to everyone with the best tools and training available to hospital staff, caretakers and patients.



Inside the single-use kit:

- ambIT* programmable intermittent bolus pump
- Cassette with filter
- Button guard
- Carrying pouch
- AA batteries
- Setup instructions
- Patient manual
- Clinical manual

AMBIT*+ PROGRAM

Combined with the ambIT* ambulatory infusion pump, the ambIT*+ program offers the best overall value to hospitals, staff and patients. We bundle everything you need for one price: a loaned **ambIT* programmable intermittent bolus pump**, day-two calling service, our extensive support and more.

How it works

Hospital staff give the pumps to patients, and patients return them to Avanos after use. We then replenish the hospital supply with the single-use kits and loaned pumps. It's efficient for staff, cost-effective for the hospital and great care for patients.

Day-two calling service

Ensure the best patient care with a patient engagement platform. Our day-two calling service provides you with real patient feedback data and can be customized to meet your patients' needs.

What are the benefits?

- No capital purchase
- A cost-effective option compared to other disposable, single-use pumps



AMBIT* PROGRAMMABLE INTERMITTENT BOLUS PUMPS

Meet patients' needs by combining basal, patient-controlled bolus, regularly occurring doses and manual bolus options to achieve a patient-specific dose of targeted pain management that matches the therapy prescribed.

Why choose ambIT*?

- Maintain a volumetric accuracy of ± 6 percent
- Connect to any IV bag up to 1000 mL of local anesthetic
- Adjust modes to Patient-Controlled Bolus, PIB or both
- Set basal rate between 0.1-20 mL/hr
- Set bolus between 0.1-50 mL/lockout time
- Set intermittent bolus rate up to 50 mL/dose time
- Provide education, training and 24/7 support to your staff

Note: A bolus or intermittent bolus is delivered at a rate of 210 mL/hr

Ensure patient safety and effectiveness with ambIT*:

- The use of the ambIT* Pain Control System, as part of a multimodal approach for pain management, can eliminate the need for opioids.
- ambIT* Electronic Pain Control Pumps allow physicians to tailor solutions to a specific patient need.
- The ambIT* Continuous pump supports the administration of high-risk medication in the comfort of the patient's home, while ensuring safety and effectiveness.
- The simple choice for customized opioid-sparing protocols: ON-Q* + ambIT.

1-800-444-2728 | ambitpump.com